

Bonhawa Heated Breathing Set



The Heated Breathing Circuit Set helps provide **continuous humidity** to dry gas and **maintains heat** and **humidification** of the gas during high flow oxygen therapy.

Breathing Circuit Components



① Water Chamber



② Water Chamber Adapter



③ Heated Breathing Circuit

Catalog No.: HBCK01

The Heated Breathing Circuit Set contains a single use intelligent water chamber adapter, a water chamber and a heated breathing circuit.

Water Chamber



| | |
|-------------------------------|------------------------|
| Duration of use | 7 Days, Single-Patient |
| Patient type | Pediatric/Adult |
| Fill Water Type | Automatically |
| Connection | Proprietary |
| Maximum Water Capacity | 150ml |
| Maximum Peak Flow | 180ml/min |
| Compatibility | ISO 80601-2-74 |

Heated Breathing Circuit



| | |
|------------------------|--|
| Duration of use | 7 Days, Single-Patient |
| Patient type | Pediatric/Adult |
| Circuit Length | 200 ± 20 cm |
| Compatibility | ISO 5367, ISO 5356, ISO 80601-2-74, IEC 60601-1, IEC 60601-1-2 |

Contact Us

Telesair, Inc.

199 Technology Drive, Suite 110

Irvine, CA 92618 USA

W: www.telesair.com

E: info@telesair.com

T: +1 949.570.3553